CRANIAL COVED Cranial Burr Hole Cover

Avoid **Postoperative Skin Concavities** Preserve Patient Cosmesis

ULTRA LOW PROFILE Excellent cosmetic results VARIED INDICATIONS Stand-alone or craniotomy burr holes FAST & EASY Sterile and instrument-free




STAND-ALONE OR COMBINED IN A CRANIOTOMY

Designed to cover burr holes in diverse cranial surgical procedures.

Biopsies · Neuroendoscopies · Craniotomies







ULTRA LOW PROFILE AND EXCELLENCE IN COSMESIS

Very smooth profile (<0.3 mm) and high adaptability.



Pterional area



Frontal burr hole covered

PEEK-OPTIMA® POLYMER

Reduced risk of soft tissue damage: smart design and less aggressive material.

......

CT/MRI ARTIFACT-FREE


Artifact-free in X-Ray, CT scan and MRI.

FAST & EASY APPLICATION

Supplied sterile and ready to use. Completely instrument-free: no screws needed.

1. PRODUCT INTRODUCTION **PULL & TIGHTEN 2.** REMOVE THE NON-IMPLANTABLE PARTS **3.** SELF CUTTING

www.neosurgery.com



AVAILABLE IN 2 SIZES

Adapts to surgeons' needs.

NEOS SURGERY S.L. (BARCELONA)- SPAIN

MANUFACTURER



FDA 510(k) cleared